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Update within deadline 1 month

Mandatory Yes No

Affected Warranty Contract Others

Update by Customer Informatics RSC Apps CSE

Remote update Yes No

Intranet download available Yes No

Update material required Yes No

Material free of charge Yes No No, credit if returned

Return of parts Yes No

Application training needed Yes No Recommended

Application training time 0 h

Estimated completion time 1/4 h Number of CSEs: 1

Remarks: N/A

Scope

Material number	11254336, 11284381, 11284500
Software version	NX2: Below VA10C, VA11C NX3: Below VA10G, VB20C
Component material number	N/A
Related to Update Instructions	N/A

Change reference no.: 690283
 Name: Yang, Jaehyuk
 Department: SHS CS US SLM

Document Version

Siemens Healthcare GmbH reserves the right to change its products and services at any time.

In addition, manuals are subject to change without notice. The hardcopy documents correspond to the version at the time of system delivery and/or printout. Versions to hardcopy documentation are not automatically distributed.

Please contact your local Siemens Healthineers office to order a current version or refer to our website <http://www.siemens-healthineers.com>.

Disclaimer

Siemens Healthcare GmbH provides this documentation "as is" without the assumption of any liability under any theory of law.

The content described herein shall be used by qualified personnel who are employed by Siemens Healthcare GmbH or one of its affiliates or who are otherwise authorized by Siemens Healthcare GmbH or its affiliates to use such documents.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens Healthcare GmbH or one of its affiliates are not entitled to use this documentation without prior written authority.

Copyright

"© Siemens Healthineers, 2019" refers to the copyright of a Siemens Healthineers entity such as:

Siemens Healthcare GmbH - Germany

Siemens Shenzhen Magnetic Resonance Ltd. - China

Siemens Shanghai Medical Equipment Ltd. - China

Siemens Healthcare Private Ltd. - India

Siemens Medical Solutions USA Inc. - USA

Siemens Healthcare Diagnostics Inc. - USA and/or

Siemens Healthcare Diagnostics Products GmbH - Germany

1	Preparation	4
1.1	General Information	4
1.1.1	Systems/Products Affected	4
1.1.2	Reason for the Update	4
1.1.3	Prerequisites	5
1.1.4	Special Tools	5
1.2	Material Information	6
1.2.1	Ordering Information	6
1.2.2	Content of the Update Material	6
1.2.3	Return of Parts	6
2	Update	7
2.1	Work Steps	7
2.1.1	Final Check	7
2.1.2	Final Work Steps	7
2.1.3	Customer Information	7
2.1.4	Changes to Previous Version	7
3	Completion Protocol US038/18/S	8
3.1	ACUSON NX Series CSAN Update	8
3.1.1	Reason for the Update	8
3.1.2	Protocol	9
4	List of Hazard IDs	10

1.1 General Information

1.1.1 Systems/Products Affected

ACUSON NX2 and NX2 Elite systems with system software versions VA10A, VA10B, VA11A and VA11B.

ACUSON NX3 and NX3 Elite systems with system software versions VA10A, VA10B, VA10C, VA10D, VA10E, VA10F, VB20A and VB20B.

1.1.2 Reason for the Update

When does this issue arise and what is the potential risk to health?

As part of our ongoing quality initiatives, internal testing revealed a situation where under specific imaging conditions during pulsed wave Doppler, there is a possibility for transducers to exceed the acoustic output power (AOP) values defined in Track 3 of the FDA Guidance for Industry and Staff: Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers and the Output Display Standard in IEC 60601-2-37.

The probability of occurrence of this issue is remote. During routine clinical use, input values to the calculations for acoustic output power are updated each time the user presses the UPDATE key for real-time acquisition. If the issue were to occur, there is an extremely remote likelihood these imaging conditions could result in a burn or tissue damage from cavitation. We have no reports of injuries or adverse events due to this issue.

What steps can the user take to avoid potential risks associated with this issue?

To avoid a situation where incorrect acoustic output power values are used during pulsed wave Doppler, always press the **UPDATE** key prior to acquiring a real-time (live), pulsed wave Doppler spectrum. The ultrasound system updates the acoustic output power values to the correct values after each press of the key.

Also, do not use a simultaneous format for continuous acquisition. You must press the **UPDATE** key prior to acquiring a real-time (live), pulsed wave Doppler spectrum. For additional information on simultaneous format, refer to Chapter A4 in your Features and Applications Reference manual.

Observe the following instructions provided for acoustic output power in Chapter 2 of your Instructions for Use.

- Acoustic Output - Mechanical and Thermal Indices

Warning: Ultrasound procedures should be used for valid reasons, for the shortest period of time, and at the lowest mechanical/thermal index setting necessary to produce clinically acceptable images.

The ultrasound system incorporates an output display of Mechanical and Thermal Indices to allow you to monitor, and to limit, the amount of ultrasound energy that is transferred to the patient.

Note: For systems distributed in the United States of America, refer to the Medical Ultrasound Safety ultrasound education program brochure produced by the AIUM that is shipped with the ultrasound system.

How will this issue be resolved?

A software update which will address this issue is currently being developed. When this software becomes available, a Siemens Healthineers customer service representative will contact you to schedule your system update. If you have any further questions, please contact your Customer Care Center.

Please share this information with all personnel within your organization who need to be aware of this issue.

1.1.3 Prerequisites

The latest version of this document can be downloaded from the (→ Services Knowledge Base (SKB), <https://skb.healthcare.siemens.com>).

1. Use **PKI card** or **Windows Login** to login.
2. Check the document type "**UI**".
3. Type in the update number **US038/18/S** in the search field.
4. Click **Search**.

1.1.4 Special Tools

N/A

1.2 Material Information

1.2.1 Ordering Information

There is no material kit for this Update program.

The document below is included with this UI in the SKB.

Quantity.	Update material number	Description
1	11507581	NX2/NX3 Customer Safety Advisory

1.2.2 Content of the Update Material

N/A

1.2.3 Return of Parts

N/A

2.1 Work Steps

To the Regional Director Customer Service

- Review the complete Update Instruction.
- Identify affected systems within your region, inform BA US if any affected systems are not already included in this update.
- Distribute the attached user's country specific language Advisory Notification PN 11507581 in form of a certified mail. The reception needs to be acknowledged. Keep a record within the country to prove the letter was sent & received.

2.1.1 Final Check

N/A

2.1.2 Final Work Steps

1. If necessary, fill out and make a copy of the attached "Completion Protocol" and file it in the corresponding binder.
2. The update shall be reported by authorized personnel, either via an intranet application or by automated data transfer (only for countries connected via update handling interface).

The update is not considered to have been performed until the IVK structure in the Installed Base is updated.

2.1.3 Customer Information

N/A

2.1.4 Changes to Previous Version

N/A

3.1 ACUSON NX Series CSAN Update

3.1.1 Reason for the Update

When does this issue arise and what is the potential risk to health?

As part of our ongoing quality initiatives, internal testing revealed a situation where under specific imaging conditions during pulsed wave Doppler, there is a possibility for transducers to exceed the acoustic output power (AOP) values defined in Track 3 of the FDA Guidance for Industry and Staff: Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers and the Output Display Standard in IEC 60601-2-37.

The probability of occurrence of this issue is remote. During routine clinical use, input values to the calculations for acoustic output power are updated each time the user presses the UPDATE key for real-time acquisition. If the issue were to occur, there is an extremely remote likelihood these imaging conditions could result in a burn or tissue damage from cavitation. We have no reports of injuries or adverse events due to this issue.

What steps can the user take to avoid potential risks associated with this issue?

To avoid a situation where incorrect acoustic output power values are used during pulsed wave Doppler, always press the **UPDATE** key prior to acquiring a real-time (live), pulsed wave Doppler spectrum. The ultrasound system updates the acoustic output power values to the correct values after each press of the key.

Also, do not use a simultaneous format for continuous acquisition. You must press the **UPDATE** key prior to acquiring a real-time (live), pulsed wave Doppler spectrum. For additional information on simultaneous format, refer to Chapter A4 in your Features and Applications Reference manual.

Observe the following instructions provided for acoustic output power in Chapter 2 of your Instructions for Use.

- Acoustic Output - Mechanical and Thermal Indices

Warning: Ultrasound procedures should be used for valid reasons, for the shortest period of time, and at the lowest mechanical/thermal index setting necessary to produce clinically acceptable images.

The ultrasound system incorporates an output display of Mechanical and Thermal Indices to allow you to monitor, and to limit, the amount of ultrasound energy that is transferred to the patient.

Note: For systems distributed in the United States of America, refer to the Medical Ultrasound Safety ultrasound education program brochure produced by the AIUM that is shipped with the ultrasound system.

How will this issue be resolved?

A software update which will address this issue is currently being developed. When this software becomes available, a Siemens Healthineers customer service representative will contact you to schedule your system update. If you have any further questions, please contact your Customer Care Center.

Please share this information with all personnel within your organization who need to be aware of this issue.

3.1.2 Protocol

This update has been completed successfully. The customer was informed of the solution to the "Customer (Safety) Advisory Notice": **US038/18/S**.

Customer:

Customer No.: Func. Location

Material No. Serial number
of the system: of the system:

Notification:

Remark:

Remark:

Country: Site:

Performed by: Telephone:

Date: Signature:

File the protocol in the corresponding binder.

There are no Hazard IDs in this document.

- Restricted - All documents may only be used by authorized personnel for rendering services on Siemens Healthcare Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hardcopies is prohibited. Offenders will be liable for damages. All other rights are reserved.

siemens-healthineers.com

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Telephone: +49 9131 84-0
siemens-healthineers.com

Print No.: US34-105.896.01.01.02 | Replaces: n.a.
Doc. Gen. Date: 01.19 | Language: English
© Siemens Healthineers, 2019